

BESTMAN



Feeding Pump



BESTMAN

BFP-10

Feature

Microcomputer control, accurate and reliable.

Friendly interactive interface , progressive guidance for procedures and operations.

Automatic alarm , safe and reliable

Tubing can quickly connect with bottled or bagged sets.

Immediate memory function to ensure the infusion dose, individualized speed and accuracy.

A variety of sensor technology that can detect the state of the upstream and downstream flows.

Automatic perfusion function.

AC and DC power supply

Handle and sucker design ,easy to carry and fixed liquid nutrients pre-heating(for option)



Specification

Flow speed

Set range 1-400ml/hr Incremental: 1ml/hr

flow speed accuracy 10% of the selected flow rate or

0.5ml/hr, whichever is greater for sure

Infusion of the total range 1-3000ml

every time 1-1000ml

Count 1-99 times



Interval 1-20 hours

Rinse 10-50ml

Interval 1-25 hours

blocked pressure Max 15psi(103kpa)





CE0086



ISO9001



ISO13485



SHENZHEN BESTMAN INSTRUMENT CO.,LTD

www.szbestman.cn

Add: 8/F,YiFang Bldg,No.315 Shuangming Ave,Dongzhou Community, Guangming St.,Guangming Dist., Shenzhen, Guangdong,518107, China

Tel: +86-755-26411423

Fax: +86-755-26495167